



Orbital rotary shaker:

Speed range: 30-300 rpm

Amplitude of rotation: 20~25 mm

Permissible ambient temperature: +5 ~ 50 C

Permissible relative moisture: N/A

Timer & alarm: 96 hours

Display: LCD

Flask holder: Universal

Holder: for 1000 ml Max

Platform size: 420mm x 420 mm

Outer Dimension (w x d x h): 600 mm x 600mm x 600mm

Power supply: AC 230V

Price: BD Tk.185,000.00 (5~ 50 C)

Price: BD Tk.145,000.00 (30~50 C)